












<b>AVOCADO</b>	<b><u>OBSERVATION REPORT</u></b>	
	Temperature : 28°C	Humidity : 80%
	Report Prep. Date	22.07.2021

Sr. No	Observation	Control	Treated with Sachets only	Treated in KIF/MAP Bags + Sachets
1	Snap			
	Colour	Dark Green	Dark Green	Dark Green
	Stage	Unripe	Unripe	Unripe
	TSS	-	-	-
	Hardness (FPT)	Hard	Hard	Hard
	Defects	None	None	None




## After 24 Hours

2	Snap			
	Colour	Greenish Black	Dark Green	Dark Green
	Stage	Semi Ripened	Unripe	unripe
	TSS	-	-	-
	Hardness (FPT)	Soft as compared to treated Samples	Hard	Hard
	Defects	None	None	None




## After 48 Hours

<b>3</b>	<b>Snap</b>			
	<b>Colour</b>	<b>Black</b>	<b>Greenish Black</b>	<b>Greenish Black</b>
	<b>Stage</b>	<b>Semi Ripened</b>	<b>Semi Ripened (3 out of 4 pcs)</b>	<b>Semi Ripened (2 out of 4 pcs)</b>
	<b>TSS</b>	-	-	-
	<b>Hardness (FPT)</b>	<b>Soft</b>	<b>Soft (3 out of 4 pcs)</b>	<b>Soft (2 out of 4 pcs)</b>
	<b>Defects</b>	<b>None</b>	<b>None</b>	<b>None</b>




## After 72 Hours

<b>4</b>	<b>Snap</b>			
	<b>Colour</b>	<b>Black</b>	<b>Black</b>	<b>Greenish Black</b>
	<b>Stage</b>	<b>Ripened</b>	<b>Ripened</b>	<b>Semi Ripened (2 out of 4 pcs)</b>
	<b>TSS</b>	-	-	-
	<b>Hardness (FPT)</b>	<b>Soft</b>	<b>Soft</b>	<b>Soft (2 out of 4 pcs)</b>
	<b>Defects</b>	<b>None</b>	<b>None</b>	<b>None</b>

## After 96 Hours

<b>5</b>	<b>Snap</b>			
	<b>Colour</b>	<b>Black</b>	<b>Black</b>	<b>Black</b>
	<b>Stage</b>	<b>Ripened</b>	<b>Ripened</b>	<b>Ripened (3 out of 4 pcs)</b>
	<b>TSS</b>	-	-	-
	<b>Hardness (FPT)</b>	<b>Soft</b>	<b>Soft</b>	<b>Soft (3 out of 4 pcs)</b>
	<b>Defects</b>	<b>None</b>	<b>None</b>	<b>None</b>

## After 120 Hours

<b>6</b>	<b>Snap</b>			
	<b>Colour</b>	<b>Black</b>	<b>Black</b>	<b>Black</b>
	<b>Stage</b>	<b>Ripened</b>	<b>Ripened</b>	<b>Ripened (3 out of 4 pcs)</b>
	<b>TSS</b>	-	-	-
	<b>Hardness (FPT)</b>	<b>Soft</b>	<b>Soft</b>	<b>Soft (3 out of 4 pcs)</b>
	<b>Defects</b>	<b>None</b>	<b>None</b>	<b>None</b>

Keep It Fresh LLP - 18 Sec 6 Manesar, Gurugram, India

USA - 14310 Gannet St , La Mirada, California

[info@keep-it-fresh.com](mailto:info@keep-it-fresh.com) [www.keep-it-fresh.com](http://www.keep-it-fresh.com)

USA | Italy | France | India | China | Mexico | Germany

